PRODUCTINFORMATIE GENU COMFORT

Dutch

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PRODUCT INFORMATION GENU COMFORT

English

Description: The Genu Comfort is a soft knit knee bandage to restore mobility, relieve pain, promote healing and quickly resolve swelling. Protects the joint without restricting mobility.

Sizing:

MAAT	KNEE SIZE	ART.NR.
S	31 - 34 cm	PO 31701 S
M	34 - 37 cm	PO 31701 M
L	37 - 40 cm	PO 31701 L
XL	40 - 44 cm	PO 31701 XL
XXL	44 - 48 cm	PO 31701 XXL

Indications: Mild knee instability, swelling and pain

Manual:

- Push the Genu Comfort over your knee and make sure that the gray padding is on the front of the kneecap.
- Make sure the clamp is not too tight. Excessive compression can lead to decreased blood flow.

Important information and warnings:

This medical product is intended for use by one individual patient. The product should be used in accordance with the advice and instructions of your physician. The product should be installed by a trained professional the first time.

In cases where the product is recommended by the doctor to be used at night, it must be ensured that it does not interfere with blood circulation. The product should be detached in case of numbness and removed if necessary. Consult your doctor in case of persistent complaints.

The product should not be exposed to excessive heat or flame. Contact with oily or acidic substances or lotions containing pomades should be avoided. This product must not be used outside its scope or combined with any other product. It should only be used according to the original indications.

Daily use and duration of application should be determined by the physician. Excessive and prolonged use can cause muscle weakness. Always consult your doctor for the correct wearing schedule.

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